


EMG recording and analysis report

Generated on: 2026-09-06 07:06

Test identifier: P01

File: ejemplo.edf



Metrics

Metric	Value
Duration	23.3 s
Analysed window	116.70-140.00 s of 140.0 s
Global RMS	0.0876 mV
Mean frequency (MNF)	142.3 Hz
Median frequency (MDF)	126.5 Hz
iEMG	0.806 mV·s
Fatigue evidence	Not conclusive — the trend does not fit (+0.20 Hz/s, $R^2=0.00$). Fatigue needs a con

Calibration (maximal voluntary contraction)


Muscle	Reference	Source
FCR	0.218 mV	calibration in this recording
ECR	0.173 mV	calibration in this recording

Calibration repetitions

Muscle	Repetition	Value
FCR	1	0.218 mV
FCR	2	0.179 mV
FCR	3	0.160 mV

FCR	4	0.140 mV
FCR	5	0.130 mV
FCR	6	0.148 mV
ECR	1	0.120 mV
ECR	2	0.120 mV
ECR	3	0.173 mV
ECR	4	0.131 mV
ECR	5	0.136 mV
ECR	6	0.119 mV

Contractions



#	Start (s)	Duration (s)
1	0.5	0.42
2	1.1	1.04
3	2.4	0.99
4	3.7	1.01
5	4.9	1.12
6	6.2	1.75
7	8.6	0.84
8	10.2	0.50
9	11.0	1.07

10	12.1	0.47
11	12.8	0.54
12	14.1	0.77
13	15.4	0.97
14	16.7	0.78
15	18.1	0.62
16	18.7	2.55
17	21.4	0.95

Co-activation (Falconer-Winter)

Whole recording — accept the fragments for one value per window

Window	FCR — Mean activation (%)	EMG — Mean activation (%)	CMV — Co-activation index
Whole recording (0.0–23.3 s)	14	11	30 %

Configuration used

Parameter	Value
Sampling rate	1000 Hz
Channel	FCR
Band-pass	20.0-450.0 Hz
Notch (mains)	50.0 Hz
Envelope (low-pass)	5.0 Hz
RMS window	50.0 ms
Device	not stored in the EDF